



Nidra is a drug-free, wearable device that uses novel TOMAC technology to stimulate specific fibers of the peroneal nerve, reducing RLS symptoms and improving sleep quality.



Compact and easy-to-use for **fast, long-lasting relief**^{1,2,3} whenever you need it

 **Personalized** intensity levels to treat your symptoms

 **Easy-to-read** LED treatment level indicators

 **Simple to control** intensity and duration

 **High-capacity,** rechargeable battery

 High-performance, flexible fabric for a **comfortable** fit

 Charge Dispersion Interface (CDI) patches **maximize comfort and safety**

Finally, a **breakthrough** in RLS therapy

Using TOMAC technology, **Nidra** helps alleviate even severe RLS symptoms.^{1,2}

How does it work?

Think about the relief you get from walking or moving your legs. Nidra is similar, but the best part, your legs don't actually move.

Nidra sends signals through nerves and muscles in your legs to activate the same pathways involved with leg movements.

This means you experience comfortable relief that can be delivered before or during sleep.



Getting started with **Nidra** is easy!



Personalize
the therapy levels
during your Personalized
Activation session.




Activate
the bands at
the onset of RLS
symptoms.



Experience relief
in less than
30 minutes³

Scan to learn more



INDICATION

The Nidra Tonic Motor Activation (TOMAC) System is intended to reduce symptoms of primary moderate-severe Restless Legs Syndrome and to improve sleep quality in adults refractory to medications. **Rx Only.** Caution: Federal law restricts this device to sale by or on the order of a medical professional.

Please visit www.nidraRLS.com/isi for complete Safety information including all contraindications and potential complications.